


No.	PN	Desc.	UoM	QTY.
1		Rod,3/8", Nylon,White	Ea	1
2		.40 thru holes	Ea	2
3		1/4"x3/8"x2" EVA tubing	Ea	3
4		50L Bag,RCFM-2001A	Ea	1
5		Adapter 1/4"-3/8"	Ea	2
6		TPE Tubing 3/8"x5/8"x36"	Ea	2
7		Large Pinch Clamp	Ea	2
8		MPC,Female,3/8"	Ea	2
9		MPC, Male Plug	Ea	2
10		Adapter,1/8"-1/4"	Ea	1
11		TPE tubing, 1/8"x1/4"x24"	Ea	1
12		Luer,Female,1/8"	Ea	1
13		Luer, Cap	Ea	1
14		Small Pinch Clamp	Ea	1



Quote Drawing

Dimensional Tolerances
 All tubing dimension as specified.
 Fractional dimensions (+/- 1/16")
 0.0 (+/-0.1"), 0.00 (+/-0.05), 0.000(+/-0.005)

Optimum Processing Inc.
 336 Bon Air Center #128
 Greenbrae, CA 94904
 ph: 415.722.5747

50L 36xMPC/24xLuer Container Assy.

SIZE	SCALE	Third angle projection
A	NTS	

Drw: DM Smith Date 4/26/09
 This specification contains proprietary information. No reproduction of this document is allowed without approval or consent of Optimum Processing Inc.

50L Bioprocess Flexible Container, 3 Port, 2x36" MPC, 1x24" Luer TPE lines

Rev
2